

COVER ALL YOUR BASES

GermiCide³

DIN 02317842



Dental & Medical, Equipment *Disinfectant Cleaner*

- 3 Minute broad spectrum kill
- Disinfectant + cleaner
- Multi-surface
- Bactericidal, Fungicidal, Virucidal, Tuberculocidal

Cover all your surfaces

Dental • Medical • Janitorial • Sanitation
EMS • Veterinary • Pet Care Facilities
Body Art Studio • Athlete Centres

Product		GermiCide ³	Optim 33 TB	CaviWipes	Sani Cloth Plus	TB Minuteman
Broad Spectrum	Bactericidal	✓	✓	✓	✓	✓
	Virucidal	✓	✓	✗	✗	✓
	Fungicidal	✓	✗	✓	✗	✓
	Tuberculocidal	✓	✓	✓	✓	✓
Surface Compliance	Plastic	✓	✓	✓	✓	✗
	Smooth/Textured Vinyl	✓	✓	✓	✗	✗
	Copper/Brass	✓	✗	✓	✓	✓
	Stainless Steel	✓	✓	✓	✓	✓
	Painted Surfaces	✓	✓	✓	✗	✓
	Plexiglas	✓	✓	✓	✗	✗
	Glass	✓	✓	✓	✓	✗
	Brushed Aluminium	✓	✗	✓	✗	✓
Safety	Biodegradable	✓	✓	✓	✓	✓
	Refillable Pouch	✓	✗	✗	✗	✓
	Health Rating "0"	✓	✓	2	1	N/A
	Flammability "0"	✓	✓	3	3	2
	Reactivity Rating "0"	✓	✓	✓	✓	✓

For more information on GermiCide³ or to order please contact your District Sales Manager or call Germiphene at 1.800.265.9931. Visit us at www.germiphene.com.



GermiCide³

Effectively kills the following microorganisms at room temperature: Mycobacterium bovis BCG, Duck Hepatitis B Virus (DHBV), Hepatitis C Virus (HCV), Rhinovirus ATCC VR-286, Poliovirus 1 ATCC VR-58, Herpes Simplex II Virus (HSV 2) G Strain, Human Immunodeficiency Virus (HIV-1), Influenza Virus Strain A2/Hong Kong, Staphylococcus aureus ATCC 6538, Staphylococcus aureus (MRSA), Enterococcus faecium (VRE), Pseudomonas aeruginosa ATCC 15442, Salmonella enterica ATCC 10708, Trichophyton mentagrophytes ATCC 9533.

Use on both hard and soft surfaces including: plastic, ceramic, stainless steel, chrome, acrylic, leather, vinyl, glass, Plexiglas®, and baked-on surfaces.

Compatible with most impression materials and dentures.

Non-sticky and leaves no hazardous chemical residue when used as directed.

Non-flammable. No special shipping or storage is required.


Safely dispose via standard waste disposal methods.

Active Ingredients:

- n-Alkyl dimethyl benzyl ammonium chloride.....0.154%
- n-Alkyl dimethyl ethylbenzyl ammonium chloride.....0.154%
- Isopropanol.....21.000%

Classes of Microorganisms* | GermiCide³ Disinfection

HARDEST to kill



Bacteria with spores (*Bacillus subtilis*, *Clostridium tetani*, *C. difficile*, *C. botulinum*)

Protozoa with cysts (*Giardia lamblia*, *Cryptosporidium parvum*)

Mycobacteria *Mycobacterium tuberculosis*, *M. avium-intracellulare*, *M. chelonae*)


Non-enveloped viruses (*Coxsackieviruses*, *polioviruses*, *rhinoviruses*, *rotaviruses*, *Norwalk virus*, *hepatitis A virus*)

Fungi *Candida species*, *Cryptococcus species*, *Aspergillus species*, *Dermatophytes*)

Vegetative bacteria (*Staphylococcus aureus*, *Salmonella typhi*, *Pseudomonas aeruginosa*, *coliforms*)

Enveloped viruses (*Herpes simplex*, *varicella-zoster virus*, *cytomegalovirus*, *Epstein-Barr virus*, *measles virus*, *mumps virus*, *rubella virus*, *influenza virus*, *respiratory syncytial virus*, *hepatitis B and C viruses*, *hantaviruses*, and *human immunodeficiency virus*)

Try both liquid + wipes to ensure all your infection control needs are met.



Lemon

Unscented

For more information on GermiCide³ or to order please contact your District Sales Manager or call Germiphene at 1.800.265.9931. Visit us at www.germiphene.com.

* Chart Source: Infection Control Guidelines: Hand Washing, Cleaning, Disinfection and Sterilization in Health Care (CCDR 2458)

