

Comparing Intelect RPW Systems



INTELECT RPW	INTELECT RPW 2
Resistive touch screen Slower processing, need recalibration, lower image quality	Capacitive touch screen Very reactive, no calibration, higher image quality
D-Actor hand piece : pressure range 1-5 BARS	Falcon hand piece: Shorter and lighter better comfort. Pressure range 0.3-5 BARS
Only D Actor transmitters	Spine and Peri Actor transmitters for spine and fascial treatments
V Actor up to 35Hz	V Actor HF up to 50Hz (better to reduce muscle tone)
Compressor 1 st generation Larger in size, heavy, no condenser, louder	Compressor latest generation more compact, less weight, less nosierreduced noise.
No possibility to export treatment data	Blue tooth export and import data (Not available at launch, to come with Software V2)
No condenser: risk of humidity accumulation in the hand piece	New condenser implemented, removing humidity from the hand piece for longer life
Space size: occupied surface 0,35 m ²	Space size: occupied surface 0,17 m ² · Half the surface

