

## TENS COMPARAISON CHART

Description	BioMed 2000 XL	BioMed 2000 XT	Impulse TENS D5	BioStim® Plus	BioStim® LX	BioStim® M7
Modes of Operation	3	3	6	4	4	7
Amplitude (Milli Amps)	1 to 98	1 to 98	1 to 98	1 to 98	1 to 98	1 to 98
Pulse Rate (Hz)	2 to 150	2 to 150	1 to 120	1 to 150	1 to 150	1 to 200
Pulse Width (µs)	50 to 250	50 to 250	25 to 250	10 to 250	10 to 250	10 to 250
Constant	YES	YES	YES	YES	YES	YES
Burst Programmed	YES	YES	YES	YES	YES	NO
Burst Programmable	NO	NO	YES	YES	YES	YES
Width Modulation	YES	YES	YES	YES	YES	YES
Rate Modulation	NO	NO	YES	NO	NO	YES
Rate / Width Modulation	NO	NO	NO	NO	NO	YES
SD Strength Duration	NO	NO	NO	NO	NO	YES
SD 2 Strength Duration	NO	NO	NO	NO	NO	YES
Battery	9 Volt (1)	9 Volt (1)	AA (2)	AA (2)	AA (4)	AA (2)
Timer	NO	YES	YES	YES	YES	YES
Low Battery Indicator	NO	NO	YES	YES	YES	YES
Digital	NO	NO	YES	YES	YES	YES

Distributed by:



**MTM**  
MÉDICAL TRONIK

190 Boul. St-Elzéar O., Laval (Qc) H7L 3N3

Tel : **450 669-8985** □ **1 800 361-0877**

Fax : **450 669-9532** □ **1 888 935-7001**

[www.medicaltronik.ca](http://www.medicaltronik.ca)