

maxill

MYCOBACTERICIDAL • VIRUCIDAL† • BACTERICIDAL • FUNGICIDAL



Minuteman NEX GEN DISINFECTANT & SPRAY



tb Minuteman NEX GEN Disinfectant is a new addition to maxill's line-up of hard surface disinfectants. When used as directed, tb Minuteman NEX GEN Disinfectant will disinfect in 1- and 3-minutes mycobacteria, small non-enveloped viruses, fungi and an array of pathogenic bacteria. In addition to disinfecting against HIV-1, tb Minuteman NEX GEN Disinfectant has the added benefit of disinfecting against Hepatitis B (DHBV) and Hepatitis C (BVDV) all in 1-minute, to assure significant reduction of risk of cross-infection from these three commonly feared bloodborne viral pathogens in healthcare facilities and emergency vehicles.

The formulation for tb Minuteman NEX GEN Disinfectant has been designed using four active ingredients and supplemental ingredients to work synergistically to create an unparalleled broad-spectrum disinfectant that is ideal for repeated daily use in dental offices, hospitals, nursing homes, and EMS & police vehicles. tb Minuteman Hard Surface Disinfectant may be used on a multitude of hard surfaces to inactivate a broad spectrum of pathogenic microorganisms.

- Broad Spectrum Bactericide, Virucide, Mycobactericide and Fungicide
- Meets CDC guidelines for intermediate level surface disinfection
- Cleaner and disinfectant in one product; save money, time and storage space
- Compatible with a wide range of hard, non-porous surfaces
- Available in a convenient bottle with sprayer and an economical refill jug

Designed with *You* in mind!

Product Availability

946 mL spray bottle	#61155
4.73 L jug	#61150

Visit us online for more information, including SDS's and pricing:

www.maxill.com/nexgen

1-minute contact time

- Human Immunodeficiency virus type 1 (HIV-1)*
- Duck Hepatitis B Virus (Surrogate for Human Hepatitis B Virus)*
- Bovine Viral Diarrhea Virus (Surrogate for Human Hepatitis C Virus)*
- Rotavirus
- Human Coronavirus
- Avian Influenza A (H3N2) Reassortant virus
- Herpes Simplex Virus Type 1
- Herpes Simplex Virus Type 2

3-minute contact time

- Mycobacterium bovis* BCG
- Pseudomonas aeruginosa*
- Staphylococcus aureus*
- Salmonella enterica*
- Vancomycin Resistant *Enterococcus faecalis* (VRE)
- Methicillin Resistant *Staphylococcus aureus* (MRSA)
- Adenovirus type 5 (ATCC VR-5)
- Feline Calicivirus (Surrogate for Human Norovirus)
- Trichophyton mentagrophytes*

Made in Canada with domestic and imported parts.

The information provided is given in good faith and believed to be accurate, but no warranty or guarantee is given or implied. Conditions of use under which the product is being applied will vary from the conditions of normal laboratory efficacy testing.

MP-REF-61150B
Rev.0: 2019-11-15

maxill inc., St. Thomas, ON N5R 6C8
www.maxill.com 1-800-268-8633